

OXYCAP™

Universal Recovery Adaptor

The OxyCap™ Universal Recovery Adaptor allows for direct oxygen delivery into the airway tube. It is suitable for use with all types of Laryngeal Masks and Tracheostomy devices. Available in two variants, with or without respiratory indicator.



- Available in two variants – with or without respiratory indicator. The respiratory indicator allows for detection of CO₂ from within the airway tube
- Cost effective – reduces the need for an oxygen mask by facilitating direct oxygen delivery into the airway tube
- Suitable for use on all types of Laryngeal Masks and Tracheostomy devices on spontaneously breathing patients
- Red colour connector – allows for quick, visual identification of OxyCap™ connected patients
- Hydrophobic filter – fitted to the respiratory indicator to reduce the risk of monitor contamination
- Lightweight device – reduces the risk of accidental displacement of the airway tube
- Individually wrapped
- MRI compatible
- Single Use
- Latex Free



Laryngeal Masks



Tracheostomy Devices

Ordering details

Product code	Product description	Units per box
OC0100	OxyCap™ with respiratory indicator	50
OC0200	OxyCap™ without respiratory indicator	50