

Stat-Check CO₂ Indicators



The STAT-Check is used as an adjunct assessment tool for verification of successful ventilation/ intubation. Ventilation, intubation and circulation are confirmed by carefully watching for colour change of the CO₂ indicator.

Every second counts!

A simple twist removes the cap and immediately activates CO₂ detection

- Both DEHP and latex free
- Can be used for <4 hours
- Ideal for emergency situations



Suitable on all patients >2.2 Kg body weight

Stat-Check CO₂ Indicators

The colour of the indicator cartridge suggests the percentage of the end-tidal CO₂ in the concentration range of 0 to 5%.

The colour of the cartridge may be compared to the label to give an approximation of end-tidal CO₂. Lack of colour change may indicate improper intubation.

Colour indication of CO₂ Presence



Active Resuscitation / Ventilation



Product code	Product description	Units per box
SC200	STAT-Check CO ₂ Indicator	10