

BactiClear[®]

ANTIMICROBIAL
Disposable Suction Liner
Range



VacSax[®]

Our **commitment** is to **quality**



Certificate No. 5482
ISO 9001:2000 & EN ISO 13485

BactiClear®

ANTIMICROBIAL

Disposable Suction Liner Range

The only system where all external surfaces are totally antimicrobial, designed to minimise the risk of cross infection.

Liners are available in 1 Litre, 2 Litre and 3 Litre versions. With canisters which fit any bracket within the hospital, suitable for ALL Clinical Areas.

Use as a completely closed system with the standard or PreGelled suction liners. Sample port versions for the addition of solidifying gel or for using with a sample sock are also available. All our liners with the exception of the cascade versions have a hydrophobic bacterial shut off filter for complete overflow protection.

Our unique patented designed liner is folded for ease of storage and begins to inflate immediately the vacuum is attached. With the integral carry handles, the liner allows for the safe removal from the canister, standing upright, whilst full of fluid, this makes

BactiClear® Antimicrobial Disposable Suction Range from **VacSax®** the product of choice



Certificate No. 542
ISO 9001:2000 & ISO 13485:2003



All VacSax products are CE marked



Class 1 Measuring Applies only to CMV
0473

Determination of antibacterial activity using modified ISO22196
Certificate Number -MSSL-LRN-1017863.1A-ISO-22196
Certificate Number -MSSL-LRN-1017863.1B-ISO-22196



Cascade Stand detail



8V 2 Litre Cascade Stand set up



Cascade Stand detail




BactiClear® Range
(see reverse for specifications)



1, 2 & 3 Litre Liners and Canisters.



Manufactured in the UK to exacting standards.



BactiClear®

ANTIMICROBIAL

Disposable Suction Liner Range

VacSax BactiClear® Antimicrobial Suction Liner with Hydrophobic Bacterial Filters

Part No.	Description	Quantity
9910-340	1 Litre Antimicrobial Liner	box of 25
VAL-102	1 Litre Antimicrobial Wide Bore Liner	box of 30
VAL-201	2 Litre Antimicrobial Liner	box of 30
VAL-301	3 Litre Antimicrobial Liner	box of 20
VAL-202	2 Litre Antimicrobial Wide Bore Liner	box of 30
VAL-104	1 Litre Antimicrobial Liner with sample port	box of 30
VAL-204	2 Litre Antimicrobial Liner with sample port	box of 30
VAL-304	3 Litre Antimicrobial Liner with sample port	box of 20



VacSax BactiClear® Antimicrobial Cascade Liners including Cascade/Tandem Tubing

Part No.	Description	Quantity
VAL-203	2 Litre Antimicrobial Cascade Liner including Cascade Tubing	box of 30
VAL-303	3 Litre Antimicrobial Cascade Liner including Cascade Tubing	box of 20
VAL-205	2 Litre Antimicrobial Sample Port Cascade Liner including Cascade Tubing	box of 30
VAL-306	3 Litre Antimicrobial Sample Port Cascade Liner including Cascade Tubing	box of 20




VacSax BactiClear® Antimicrobial PreGelled Liners with Hydrophobic Bacterial Filters

Part No.	Description	Quantity
9810-101	1 Litre PreGelled Antimicrobial Liner	box of 25
GAL-211	2 Litre PreGelled Antimicrobial Liner	box of 30
GAL-311	3 Litre PreGelled Antimicrobial Liner	box of 20
GAL-212	2 Litre Wide Bore PreGelled Antimicrobial Liner	box of 30



VacSax BactiClear® Antimicrobial PreGelled Cascade Liners including Cascade/Tandem Tubing

Part No.	Description	Quantity
GAL-213	2 Litre PreGelled Antimicrobial Cascade Liner including Cascade Tubing	box of 30
GAL-313	3 Litre PreGelled Antimicrobial Cascade Liner including Cascade Tubing	box of 20



VacSax Limited

Western Wood Way, Langan Science Park,
Plymouth, Devon PL7 5BG, U.K.

Tel: +44 (0)1752 337000

Fax: +44 (0)1752 337600

Email: sales@vac.sax.com

www.vacsax.com



Certificate No. 542
ISO 9001:2000 & EN 46001:1996



VacSax®

Our commitment is to quality