

NacSax[®]

Disposable Suction Irrigation Set



Certificate No. 542
ISO 9001:2008 & ISO 13485:2012



Our **commitment** is to **quality**



Disposable Suction Irrigation Set

- Unblocked channel with good flowing performance.
- Transparent hand piece, unobstructed visualization. Stainless steel probe with position function, easy to fix and draw out.
- Ergonomically designed hand piece can be easily operated by either hand.
- Clear channel, consistent diameter with self-flush function.
- Saline evacuation and filtration, Saline flow rate is adjustable.
- Flip designed rear port for easy access inside chamber, can be operated by one hand, quick close and reduce air leakage.

Item Code	Probe Dia.	Probe Length	Package
9750-001	ø5mm	330mm	pack of 10
9750-002	ø5mm	450mm	pack of 10

Material: PVC + ABS + Silicone + PE/PET + SUS304



Certificate No. 542
ISO 9001:2008 & ISO 13485:2012

VacSax Limited

Western Wood Way, Langage Science Park,
Plymouth, Devon PL7 5BG, U.K.

Tel: +44 (0)1752 337000

Fax: +44 (0)1752 337600

Email: sales@vac.sax.com

www.vacsax.com