

DUET

Portable Electric Suction Unit with battery back-up



S-Scort® Duet Portable Suction with 2 Litre VacSax Liner and Universal Canister Holder

The **S-Scort® Duet** is a powerful, portable AC suction unit with battery back-up. With airflow greater than 30Lpm at the end of a patient tube, the Duet will give you the suction needed for all constant suction and emergency care procedures.

The **S-Scort® Duet** can be conveniently positioned on stand-by anywhere you need and the optional retention bracket enables the Duet to be fastened to a Trolley or Wall.

The battery care system automatically charges and monitors itself. Leave it plugged in all the time, the battery cannot be overcharged. The automatic capacity check will warn you when the battery begins to lose capacity (eventually all batteries lose capacity and need to be replaced, an explanation of this can be found at www.sscor.com/maxcapacity.htm)

Sscor will replace a new Duet battery that does not operate to specifications within three years of a new Duet purchase.

4720-203
Universal Canister Holder



S-Scort® Duet Portable Suction with 1 Litre VacSax Liner and filter kit



DUET

Portable Electric Suction Unit with battery back-up

SSCOR, INC. DUET

- 9720-210 - S-Scort® Duet with Universal Holder and filter kit
 - 9720-211 - S-Scort® Duet with Filter Kit
 - 4720-202 - Retention Bracket (for wall or Trolley)
 - 4720-203 - Universal Canister Holder
 - 4720-209 - Filter Kit
 - 3833-134 - 1L BactiClear Antimicrobial Suction Canister
 - 9910-340 - 1L BactiClear Antimicrobial Liner - Box of 25
 - 9810-101 - 1L BactiClear Antimicrobial PreGelled Liner - Box of 25
- BactiClear Liners - Patent No. GB2497360

S-Scort® Specifications

- Performance** 30Lpm airflow at the end of patient tube:-
50mmHg to -550mmHg vacuum
- Dimensions** 433 x 229 x 134mm
- Weight** 3.9kg
- Charger** Internal automatic charger offers fast charge up to full and then gentle float to maintain charge. LED's indicate charging status and operation mode.
- Manufacturer** **SSCOR, INC.**, Sun Valley CA 91352 USA



All liners manufactured by VacSax are made from High Density PE and do not contain any PVC and produce only carbon dioxide and water upon 100% combustion helping to save the environment.

Every effort has been made to ensure the accuracy of information contained within this literature. However freedom from errors and omissions cannot be guaranteed. VacSax policy is one of continued product improvement. The company reserves the right to change the information contained herein without prior notice.

VacSax Limited is an ISO 13485 registered company



Our commitment is to quality

VacSax Limited

Western Wood Way, Langage Science Park,
Plymouth, Devon PL7 5BG, U.K.

Tel: +44 (0)1752 337000 **Fax:** +44 (0)1752 337600

Email: sales@vac sax.com

www.vacsax.com