

Absorbent Theatre Mat



Nordiska - Absorbent Theatre Mat

Suitable for use under the operating table

Dual-sided

Allows the mat to absorb, store and bind fluid through both sides

Superabsorbent Core

Absorbs and retains liquids in a safe and hygienic manner


Polypropylene non-woven fabric outer layer

Provides a fatigue resistant cover, while allowing fluid to pass through into the core

Single-use

Non-sterile

Latex-free



| | |
|---------------------|---------------------------|
| External Dimensions | 75cm x 36cm |
| Capacity | Up to 8 litres |
| Pack Size | 100 pieces (mats) per box |
| Product Code | 95922 |