

CLEAN MEDICAL PRODUCTS

Engineer For Health

CE 0123



CLEANEST AC Pulse Lavage MODEL W-205

The use of the CLEANEST and its effectiveness have been clinically proven and well documented in both orthopaedic surgeries and in wound management.

In orthopaedic surgeries, the CLEANEST uses powerful pulsed irrigation to thoroughly clean up the cancellous bone bed by washing away blood and debris, making possible a much improved bone to cement interface and subsequently longer survival of the prosthesis while also reducing infection rate and risk of fat embolism.

In various types of wound debridement from burn to trauma to chronic wounds, the CLEANEST effectively and efficiently removes more than 90% of the necrotic tissues, foreign contaminants and bacteria; facilitates a faster healing process while greatly reduces infection rate.

Guangzhou Clean Medical Products Manufacturing Corp.
No.9 Guangcong Rd, Conghua Development District, Guangdong, China
Tel: +86-20-87654060, 87627597 Fax: +86-20-87303882
Email: info@cleanmfg.com Web: www.cleanmfg.com

-1-

Adjustable high and low speed for orthopaedic and soft tissue debridement.

-2-

Powerful irrigation pressure for orthopaedic and gentle pressure for wounds.

-4-

Double packed and sterilized.

-3-

Reliable and consistent power.

-6-

Ergonomic design with minimal installation; features nozzle lock to secure nozzle in place during procedures.

-5-

Integrated tubings to avoid intertwining.



Large Splash Guard



1
Less Waste: The re-usable AC adaptor allows the riddance of any hazardous alkaline batteries and is an environmentally friendly option.

2
More Power: Extra powerful irrigation and non-diminishing power output.

Fan Spray (Knee) Nozzle
Shower Spray Nozzle



Femoral Canal (Hip) Nozzle



Femoral Brush Nozzle

REF No.	Product Description	Packaging Info.
Disposable Hand Piece		
W-205-02-A	Disposable Hand Piece with only Knee Nozzles	10 units / carton 10.5 Kg 47 x 36 x 45 cm
W-205-02-B	Disposable Hand Piece with both Femoral and Knee Nozzles	
Re-usable AC Adaptor		
W-205-01B	Re-usable AC Adaptor	8 units / carton 10.7 Kg 44 x 42 x 33 cm
Nozzles		
W-205-P1	Fan Spray (Knee) Nozzle	50 units / carton box 46 x 43 x 36 cm
W-205-P2	Femoral Canal (Hip) Nozzle	
W-205-P3	Femoral Brush Nozzle	
W-205-P4	Shower Spray (Spine) Nozzle	
W-205-P11	Large Splash Guard (12cm Diameter)	