

Swab Counting Bag

Surgical Sponge Counter Bags



Used to support operating theatre efficiency and safety as well as provide easy tracking of lap sponges, Fannin's Sponge Counter Bag is an effective tool for cleaning up sponges after surgery.

Benefits & Features

- Clear polyethylene pockets allow sponge visibility for visual estimate of patient blood loss
- Easy open pockets for the simple disposal of used sponges
- Excellent visibility for counting sponges
- Latex Free

Each Bag Can Hold:

- Ten 4x4 gauze sponges **OR** Five Lap Sponges



For more information or to make an enquiry, please contact us at:

Fannin Ltd.
Fannin House
South County Business Park
Leopards Town
Dublin 18
T: +353 (0)1 290 7000
E: info@fannin.eu

Fannin UK Ltd
Fernhouse - Suite 3
Daniels Industrial Estate
Stroud, GL5 3TJ
T: +44 (0)1453 874 130
E: sales@fannin.eu

Fannin NI
Unit 2, Heron View
Airport Road West
Belfast BT3 9LN
T: +44 (0)2890 735 588
E: nicoffice@fannin.eu

www.fannin.eu