

COAT RACK

**Removable pegs
for gowns**

PREVENTION TOOL

**Instruction display
window**

HAND HYGIENE

**Disinfectant solution
holder**

PERSONAL PROTECTIVE EQUIPMENT

**Storage baskets for glove
& surgical mask boxes**

WASTE MANAGEMENT

Refuse bag holder



INFECTION
PREVENTION
AND CONTROL

**Mobile
Isolation
Trolley
(UMP[®])**

Mobile Isolation Trolley - UMP®

Registered mark & models

Compact and modular solution, the **Mobile Isolation Trolley - UMP** is designed to ensure the rapid implementation of aseptic and protective isolation procedures.

Placed at the entrance or inside the patient's room, this educational tool alerts, informs and limits the risk of infectious transmission.



Ref. UMP COL C

Description

Available in bright yellow colour to identify the high risk infection control areas in the facility, the **Mobile Isolation Trolley - UMP®** is designed to store all the supplies required by hygiene protocols to reduce the risk of pathogen transmission in hospitals between patient, healthcare workers and visitors.

► **Mobile Isolation Trolley - UMP®**
fully equipped with removable accessories:

- 2 metal storage baskets, stackable (L. 380 x l. 315 x H. 100mm)
- 1 metal shelf (L. 350 x l. 280mm)
- 1 universal hydro-alcoholic solution holder (stainless steel wire)
- 1 refuse bag holder (stainless steel wire)
- 3 pegs for gown
- 1 instruction display window

Reference

Mobile Isolation Trolley - UMP®

Mobile Isolation Trolley - UMP®
Fully-equipped / Bright yellow colour

Ref. **UMP COL C**

Features

- Stainless steel frames fully welded (25 x 25mm)
Yellow Epoxy coating to avoid the risk of corrosion.

- 4 swivel castors with non marking rubber (Ø 50mm),
2 with brakes.

