




## IV Gravity Sets




- IV Gravity Sets with air vents and 15µm filters
- Available with 3-way connector and extension tubing
- Needle free access ports as optional extras
- Flow dial regulator (mL/hr) variants
- Burette Set range

## IV Gravity Sets




### 1029DF

IV Gravity Set with 15 $\mu$ m filter, air vent, 180cm tubing, rotating luer connector




### 1029DF2

IV Gravity Set with 15 $\mu$ m filter, air vent, 240cm tubing, rotating luer connector




### 1029YDF

IV Gravity Set with 15 $\mu$ m filter, air vent, 180cm tubing, Y-site, rotating luer connector



### 1029DFV1


IV Gravity Set with 15 $\mu$ m filter, air vent, 70cm tubing, back check valve connector



### 1029DFV


IV Gravity Set with 15 $\mu$ m filter, air vent, 180cm tubing, back check valve connector

## IV Gravity Sets with extension tubing




### 1029DFX

IV Gravity Set with 15 $\mu$ m filter, air vent, 180cm tubing, 3-way connector, 120cm extension tubing, rotating luer connector




### 1029DFXV

IV Gravity Set with 15 $\mu$ m filter, air vent, 180cm tubing, back check valve, 3-way connector, 120cm extension tubing, rotating luer connector



### 1029DFXMY


IV Gravity Set with 15 $\mu$ m filter, air vent, 180cm tubing, 3-way connector, needle free access port, 120cm extension tubing, rotating luer connector



### 1029DFXVMY


IV Gravity Set with 15 $\mu$ m filter, air vent, 180cm tubing, back check valve, 3-way connector, needle free access port, 120cm extension tubing, rotating luer connector

## IV Gravity Sets with Y-site needle free access ports




### 1048BMY

IV Gravity Set with 15µm filter, air vent, needle free access port, back check valve 210cm tubing, rotating luer connector



### 1048B2MY


IV Gravity Set with 15µm filter, air vent, two needle free access ports, back check valve, 260cm tubing, rotating luer connector



### 1048BMYX


IV Gravity Set with 15µm filter, air vent, needle free access port, back check valve, 260cm tubing, rotating luer connector

## IV Gravity Sets with flow dial regulator (mL/hr)



### 1017DF


IV Gravity Set with 15µm filter, air vent, 180cm tubing, flow dial regulator, rotating luer connector



### 1017DFMY


IV Gravity Set with 15µm filter, air vent, 180cm tubing, flow dial regulator, needle free access port Y-site, rotating luer connector

## Burette Sets



### 1016DFSML

Burette Set with Y site, rotating luer connector



### 1016DFMYSML

Burette Set with two needle free access ports, rotating luer connector



Needle free access ports as optional extras



Flow dial regulator (mL/hr) variants



Available with 3-way connector and extension tubing

Product code	Product description	Units per box	NHS Supply Chain code
<b>IV Gravity Sets</b>			
1029DF	IV Gravity Set with 15µm filter, air vent, 180cm tubing, rotating luer connector	50	FSB1293
1029DF2	IV Gravity Set with 15µm filter, air vent, 240cm tubing, rotating luer connector	50	FSB1413
1029YDF	IV Gravity Set with 15µm filter, air vent, 180cm tubing, Y-site, rotating luer connector	50	FSB2214
1029DFV1	IV Gravity Set with 15µm filter, air vent, 70cm tubing, back check valve connector	50	FSB1414
1029DFV	IV Gravity Set with 15µm filter, air vent, 180cm tubing, back check valve connector	50	FSB1415
<b>IV Gravity Sets with extension tubing</b>			
1029DFX	IV Gravity Set with 15µm filter, air vent, 180cm tubing, 3-way connector, 120cm extension tubing, rotating luer connector	50	FSB1339
1029DFXV	IV Gravity Set with 15µm filter, air vent, 180cm tubing, back check valve, 3-way connector, 120cm extension tubing, rotating luer connector	50	FSB1416
1029DFXMY	IV Gravity Set with 15µm filter, air vent, 180cm tubing, 3-way connector, needle free access port, 120cm extension tubing, rotating luer connector	50	FSB1417
1029DFXVMY	IV Gravity Set with 15µm filter, air vent, 180cm tubing, back check valve, 3-way connector, needle free access port, 120cm extension tubing, rotating luer connector	50	FSB1418
<b>IV Gravity Sets with Y-site needle free access ports</b>			
1048BMY	IV Gravity Set with 15µm filter, air vent, needle free access port, back check valve, 210cm tubing, rotating luer connector	50	FSB2558
1048B2MY	IV Gravity Set with 15µm filter, air vent, two needle free access ports, back check valve, 260cm tubing, rotating luer connector	50	FSB1420
1048BMYX	IV Gravity Set with 15µm filter, air vent, needle free access port, back check valve, 260cm tubing, rotating luer connector	50	FSB2109
<b>IV Gravity Sets with flow dial regulator (mL/hr)</b>			
1017DF	IV Gravity Set with 15µm filter, air vent, 180cm tubing, flow dial regulator, rotating luer connector	50	FSB1291
1017DFMY	IV Gravity Set with 15µm filter, air vent, 180cm tubing, flow dial regulator, needle free access port , Y-site, rotating luer connector	50	FSB1772
<b>Burette Sets</b>			
1016DFSML	Burette Set with Y site, rotating luer connector	20	FSB1879
1016DFMYSML	Burette Set with two needle free access ports, rotating luer connector	20	–