


# Reusable Laryngeal Mask Airways

Silicone, reinforced,  
double lumen



# Reusable Laryngeal Mask Airways

Silicone, reinforced, double lumen



- Medical grade silicone
- Validated for up to 40 uses
- DEHP free
- Latex free

Product code	Size	Patient Weight	Cuff Volume
FF30810	1	< 5 Kg	< 4ml
FF30815	1.5	5-10 Kg	< 7ml
FF30820	2	10-20 Kg	< 10ml
FF30825	2.5	20-30 Kg	< 14ml
FF30830	3	30-50 Kg	< 20ml
FF30840	4	50-70 Kg	< 30ml
FF30850	5	70-100 Kg	< 40ml