

F-31 Suction Pump

AC/DC with Lithium Battery

The New F-31, with a powerful piston motor, is a portable suction machine for Acute and Community use. It is suitable for emergency and minor surgery



1 Complete with suction regulator and vacuum gauge

3 Supplied with 1L. sterilisable jar with overflow valve

2 Flexible and suitable for a wide range of applications complete with disposable filter and tube sets

4 Long-life oil-less piston pump motor



Reliable, easy to use, quiet and safe



Can be mounted on trolley with 4 castors



Can be charged on one 12V car-plug

Technical data

Power Supply Wall:	240V - 110 50/60Hz	Max Aspiration:	- 0,80 Bar
Power Consumption:	1,5A	IP:	33
Battery Life:	100 Min.	Noise dBA:	61,2
Battery Charging Time:	90 Min.	Min. Size Length Width Height:	45x15x28 cm
Duty Time:	60 Min. On - 60 Min. Off	Max. Size Length Width Height:	45x15x30 cm
Max Aspiration Flow Direct	40 Lt/Min.	Weight:	3,2kg
Measurement (On Aspiration Vent):		Vehicle Charging Cable:	12v (Pin 2,1mm / Ø 5,5mm)

Areas for use

- Minor Surgery
- Emergency Surgery
- Hospitals and emergency vehicles
- Veterinary and GP Practices

Description

F-31 Suction Pump AC/DC with Lithium Battery

Part number

10601/vacsax kit

Vacsax Limited

Western Wood Way,
Langage Science Park, Plymouth, Devon
PL7 5BG, U.K.

Tel: +44 (0)1752 337000

Fax: +44 (0)1752 337600

Email: sales@vacsax.com

www.vacsax.com

www.fannin.eu