

# MULTI SHEATH

## DISPOSABLE MULTI-PURPOSE INSTRUMENT SHEATH



The single-use, non-sterile, multi-sheath cover is an economical and safe way to reduce contamination of reusable instruments such as a hyfrecator handpiece

### Features

- For use with a hyfrecator hand-piece and other surgical instruments
- Protects the instrument from external contamination during procedures
- Promotes speedy patient changeover
- Highly recommended for infection control
- Conveniently packed in quantities of 200 to reduce packaging waste

Code	Description	Quantity
FHCMS200	Disposable multi-purpose instrument sheath	200