



The high performer of stone extraction baskets:

StoneWrapper[®]



Perfect for removing concretions of different sizes from hard-to-reach areas of the renal pelvis system comfortably, safely and quickly.

- Perfect hybrid of grabber and stone extraction basket
- Elastic, fully articulated three-armed nitinol basket
- Easy insertion into the lower calyx and quick removal of all concretions
- Available in two sizes and with two opening options

Effective: work efficiently, save time

The innovative StoneWrapper solves the major challenges involved with flexible ureterorenoscopy: fast and safe atraumatic removal of marginal stones from the renal pelvic system.

The combination of

- a specially designed three-arm hybrid head (combination of basket and grabber),
- a new guide wire technology with high shaft stiffness and hyperflexible distal end,
- including a Tuohy-Borst adapter for easy insertion into the endoscope's working channel and as splash guard,
- a wire design designed for durability and functionality,
- an ergonomic push handle for fatigue-free working and
- excellent irrigation for best visibility

enables the reliable and fast removal of even the smallest concretions from any position.



No stone is safe

Highly technical level for best user convenience: The StoneWrapper is the ideal instrument for complex stone recovery from the most difficult positions.

The Tuohy-Borst adapter facilitates the insertion of the basket and serves as a splash protection.

Order No.	Size / FR	Length / cm	Opening Width / mm
48012553	1.5	120	8
48012554	1.5	120	11
48012555	2.0	120	8
48012556	2.0	120	11

We are urology specialists and your reliable partner when it comes to optimal user and patient comfort. Just ask and we can recommend the best product for your needs.

Manufacturer: EPflex Feinwerktechnik GmbH
D-72581 Dettingen an der Erms

Distributor:

UROTECH GmbH
Medi-Globe-Straße 1-5
D-83101 Rohrdorf OT Achenmühle
Phone: +49 (0) 8032 - 973-210
Fax: +49 (0) 8032 - 973-211
E-Mail: info@urotech.com
www.urotech.com