

Eco Sac 5mm

Retrieval through 5mm trocars

DEPLOY

- using unique introducer rod through 5mm trocar

POSITION

- black stitching to simplify identification of the Sac's mouth

RETRIEVE


- via laparoscopic incision reducing hospital stay, cost and complication




Safe and secure closure with strong monofilament drawstring













One-piece tear resistant design with unique PUR coating for maximum strength and durability




Bag can be redeployed for multiple tissue retrieval



Eco Sac 5mm Features:

-  Comprehensive solution covering laparoscopic procedures.
-  Minimally invasive; reducing length of stay, financial strain and improved patient experience.
-  Our unique coated rip-stop nylon is resistant to tearing if punctured.
-  Espiner's unique fabric designed for leak resistance.
-  Easy insertion using unique introducer rod.
-  Easy orientation via black stitching.
-  Comprehensive tissue capture with large mouth opening.
-  Safe closure with strong monofilament drawstring.
-  Includes X-ray opaque marker.
-  Up to a five-year shelf life.



Product Code	Trocar (mm)	Diameter (mm)	Length (mm)	Volume (ml)	Number Per Box	Fits through 5mm re-usable trocar	Fits through 5mm disposable trocar
EMP50ECO-R	5	38	90	50	10	✓	✓
EMP60ECO-R-10	5	70	100	60	10	✓	✓
EMP80ECO-R-10	5	70	125	80	10	✓	✓
EMP110ECO-R-10	5	55	135	100	10	✓	✓
EMP200ECO-R-10	5	70	230	200	10	✗	✓

To talk to Espiner about our cost effective Eco Sac range, please call us on +44 (0)1275 878 801 or email espiner@fannin.eu
www.espinermedical.com