



Active Cord

For connecting Medi-Globe instruments to HF generators of various manufacturers

| Catalogue No. | Description | Total Length (cm) |
|---------------|------------------------|-------------------|
| GAC-51-04-300 | ERBE 4 mm | 300 |
| GAC-51-05-300 | ERBE 5 mm | 300 |
| GAC-52-04-300 | Martin/Berchthold 4 mm | 300 |
| GAC-53-08-300 | Valleylab, 8 mm | 300 |
| GAC-54-06-300 | Olympus, 6 mm | 300 |

L-GAC-SU-01-EN / LL-0409 Vers. 1.0 (01/2021)



Accessories

Single Use

Biopsy Valves

For support of insufflation and suction

Bite Blocks

Ergonomic adaption to patients' facial outline

Cleaning Brushes

Coated brush tip

Hygienic storage and avoidance of accessory channel scratch

Active Cords



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de

Biopsy Valves & Bite Blocks – Single Use

Cleaning Brushes – Single Use



Biopsy Valve – Standard Type

Individually packed in a dispenser

| Catalogue No. | Suitable for Endoscope-Type | Color |
|---------------|--|-------|
| GAC-04-00-001 | OLYMPUS/ FUJIFILM (500 Series and newer)/Storz | Blue |
| GAC-04-01-001 | PENTAX | Red |



Biopsy Valve – with Flushing Channel

Individually packed in a dispenser

| Catalogue No. | Suitable for Endoscope-Type | Color |
|---------------|--|-------|
| GAC-04-00-002 | OLYMPUS/ FUJIFILM (500 Series and newer)/Storz | Blue |
| GAC-04-01-002 | PENTAX | Red |



Bite Block – Latex-Free

Individually packed in a dispenser box

| Catalogue No. | Description | Passage (mm) |
|---------------|--|--------------|
| GAC-01-01-001 | With tongue pressure and rubber head strap | W 19 x H 18 |



Bite Block – Latex-Free

Individually packed in a dispenser box

| Catalogue No. | Description | Passage (mm) |
|---------------|------------------------|--------------|
| GAC-01-03-001 | With rubber head strap | W 26 x H 18 |



Bite Block – Latex-Free

Individually packed in a dispenser box

| Catalogue No. | Description | Passage (mm) |
|---------------|-------------------------|--------------|
| GAC-01-03-002 | With textile head strap | W 23 x H 20 |



Bite Block – Latex-Free

Individually packed in a dispenser box

| Catalogue No. | Description | Passage (mm) |
|---------------|--------------------|--------------|
| GAC-01-01-002 | Without head strap | W 26 x H 18 |



Endoscope Cleaning Brush Set

Contains double sided channel- and port cleaning brush; non-sterile

| Catalogue No. | Description | Brush Diameter (mm) | Brush Length (mm) | Total Length (cm) | Min. Accessory Channel (mm) |
|---------------|---|---------------------|-------------------|-------------------|-----------------------------|
| GAC-03-05-230 | Brush 1: double sided Channel Cleaning Brush | 5/5 | 20/20 | 230 | 2.8 |
| | Brush 2: double sided Valve Cleaning Brush | 5/12 | 20/40 | 15 | n.a. |



Valve Cleaning Brush

Double sided; non-sterile

| Catalogue No. | Brush Diameter (mm) | Brush Length (mm) | Total Length (cm) | Min. Accessory Channel (mm) |
|---------------|---------------------|-------------------|-------------------|-----------------------------|
| GAC-03-05-012 | 5/12 | 20/40 | 15 | n.a. |



Cleaning Brush

Double sided; non-sterile

| Catalogue No. | Color | Brush Diameter (mm) | Brush Length (mm) | Total Length (cm) | Min. Accessory Channel (mm) |
|---------------|--------|---------------------|-------------------|-------------------|-----------------------------|
| GAC-02-03-230 | Green | 3/3 | 20/20 | 230 | 2.0 |
| GAC-02-05-230 | Yellow | 5/5 | 20/20 | 230 | 2.8 |
| GAC-03-05-011 | Yellow | 5/11 | 20/40 | 230 | 2.8–6.0 |

