

Hot Biopsy Forceps – Single Use



Oval
Fenestrated

Hot Biopsy Forceps

Monopolar biopsy forceps with teflon coated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GBF-11-24-180	Oval Fenestrated	2.4	180	2.8
GBF-11-24-230	Oval Fenestrated	2.4	230	2.8



L-GBF-SU-00-EN / LL-0012 Vers. 10.0 (01/2021)



Biopsy Forceps

Single Use

Sharp Standard and Serrated (Alligator) Cups

Precise and clean tissue acquisition

Versions with Teflon-Coated Metal Spiral

Friction minimized passage through accessory channel




Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de

FlexiBite Biopsy Forceps – Single Use

Cup Shapes




Bronchoscopy

Uncoated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GBF-01-18-120	Oval Fenestrated	1.8	120	2.0
GBF-02-18-120	Oval Fenestrated with Spike	1.8	120	2.0



FlexiBite Biopsy Forceps – Single Use



Gastroscopy & Colonoscopy

Uncoated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GBF-01-23-180	Oval Fenestrated	2.3	180	2.8
GBF-01-23-230	Oval Fenestrated	2.3	230	2.8
GBF-02-23-180	Oval Fenestrated with Spike	2.3	180	2.8
GBF-02-23-230	Oval Fenestrated with Spike	2.3	230	2.8



Gastroscopy & Colonoscopy

Uncoated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GBF-03-23-230	Alligator	2.3	230	2.8
GBF-04-23-230	Alligator with Spike	2.3	230	2.8



FlexiBite Biopsy Forceps – Single Use

FlexiBite Biopsy Forceps – Single Use



Bronchoscopy

Teflon-coated stainless steel spiral


Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GBF-21-18-120	Oval Fenestrated	1.8	120	2.0
GBF-22-18-120	Oval Fenestrated with Spike	1.8	120	2.0



Gastroscopy & Colonoscopy

Teflon-coated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GBF-23-23-230	Alligator	2.3	230	2.8
GBF-23-30-230	Alligator	3.0	230	3.2
GBF-24-23-230	Alligator with Spike	2.3	230	2.8
GBF-24-30-230	Alligator with Spike	3.0	230	3.2



Gastroscopy & Colonoscopy


Teflon-coated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GBF-21-18-160	Oval Fenestrated	1.8	160	2.0
GBF-21-23-180	Oval Fenestrated	2.3	180	2.8
GBF-21-23-230	Oval Fenestrated	2.3	230	2.8
GBF-21-30-230	Oval Fenestrated	3.0	230	3.2
GBF-22-18-160	Oval Fenestrated with Spike	1.8	160	2.0
GBF-22-23-180	Oval Fenestrated with Spike	2.3	180	2.8
GBF-22-23-230	Oval Fenestrated with Spike	2.3	230	2.8
GBF-22-30-230	Oval Fenestrated with Spike	3.0	230	3.2



EcoBite Biopsy Forceps – Single Use


EcoBite Biopsy Forceps – Single Use



Bronchoscopy – Biopsy Forceps Standard

Uncoated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GBF-31-18-120	Oval Fenestrated	1.8	120	2.0
GBF-32-18-120	Oval Fenestrated with Spike	1.8	120	2.0



Gastroscopy & Colonoscopy – Biopsy Forceps Standard

Uncoated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GBF-31-23-180	Oval Fenestrated	2.3	180	2.8
GBF-31-23-230	Oval Fenestrated	2.3	230	2.8
GBF-32-23-180	Oval Fenestrated with Spike	2.3	180	2.8
GBF-32-23-230	Oval Fenestrated with Spike	2.3	230	2.8



Bronchoscopy – Biopsy Forceps Standard

Teflon Coated Stainless Steel Spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GBF-41-18-120	Oval Fenestrated	1.8	120	2.0
GBF-42-18-120	Oval Fenestrated with Spike	1.8	120	2.0



Gastroscopy & Colonoscopy – Biopsy Forceps Standard

Teflon Coated Stainless Steel Spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GBF-41-23-180	Oval Fenestrated	2.3	180	2.8
GBF-41-23-230	Oval Fenestrated	2.3	230	2.8
GBF-42-23-180	Oval Fenestrated with Spike	2.3	180	2.8
GBF-42-23-230	Oval Fenestrated with Spike	2.3	230	2.8

