

Foreign Body Forceps – Single Use



Foreign Body Forceps

Teflon coated metal spiral



Alligator

Rat Tooth

Catalogue No.	Forceps Shape	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GFB-09-23-120	Rat Tooth	2.3	120	2.5
GFB-09-23-180	Rat Tooth	2.3	180	2.5
GFB-09-23-230	Rat Tooth	2.3	230	2.5
GFB-10-23-180	Alligator	2.3	180	2.5
GFB-10-23-230	Alligator	2.3	230	2.5



L-GFB-SU-00-EN / LL-0023 Vers. 6.0 (01/2021)

Polyp and Foreign Body Retrieval Devices

Single Use



Polyp and Foreign Body Retrievers

Strong & wide opening graspers
Large variety of applications

Foreign Body Forceps

Coated metal spiral
Friction minimized passage through accessory channel



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de

Polyp and Foreign Body Retrievers – Single Use

Retrieval Devices – Single Use



4-Prong Grasper

Foreign Body Retriever – Grasper with Hooks

Teflon tube

Catalogue No.	Grasper Shape	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GFB-01-25-180	4-Prong Grasper, with Hooks	2.5	180	2.8
GFB-01-25-230	4-Prong Grasper, with Hooks	2.5	230	2.8



3-Prong Grasper



4-Prong Grasper

Polyp Retriever – Grasper with Eyelets

Teflon tube

Catalogue No.	Grasper Shape	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GFB-06-25-180	3-Prong Grasper with Eyelets	2.5	180	2.8
GFB-06-25-230	3-Prong Grasper with Eyelets	2.5	230	2.8
GFB-02-25-180	4-Prong Grasper with Eyelets	2.5	180	2.8
GFB-02-25-230	4-Prong Grasper with Eyelets	2.5	230	2.8



Forked Jaws (2:2) Grasper


Foreign Body Retriever – Grasper with Forked Jaws

Teflon tube

Catalogue No.	Grasper Shape	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GFB-03-25-180	Grasper with 2:2 Jaws	2.5	180	2.8
GFB-03-25-230	Grasper with 2:2 Jaws	2.5	230	2.8



6-Wire Basket



4-Wire Basket

Foreign Body Basket

Teflon tube

Catalogue No.	Grasper Shape	Basket Diameter (mm)	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GFB-07-25-425	4-Wire Basket	25	2.5	230	2.8
GFB-07-25-625	6-Wire Basket	25	2.5	230	2.8

