






## Fitting Lithotripsy Components - Injectable Mechanical Lithotripsy

Catalogue No.	Description
GML-99-01-000	<p><b>StoneBuster Lithotripsy Handle</b></p> <p>Reusable, user-friendly twist &amp; pull technology, 1-push safety system</p> 
GML-91-26-180	<p><b>Injectable Lithotripsy Spiral</b></p> <p>Reusable, outer diameter 2.6 mm, working length 180 cm</p> 
GML-90-34-080	<p><b>Emergency Lithotripsy Spiral</b></p> <p>Reusable, outer diameter 3.4 mm, working length 80 cm</p> 

L-GML-SU-01-EN / LL-0028 Vers. 7.0 (01/2021)



# StoneBuster Lithotripsy System

Single Use

- Combination Concept**  
Components for planned and emergency lithotripsies
- Simple and Safe Handling**  
Convenient handling for component assembly
- Lithotripsy Metal Spiral with Injection Port**  
For cannulation, catching and lithotripsy without component change




**Medi-Globe GmbH**  
 Medi-Globe-Straße 1-5 • D-83101 Rohrdorf OT Achenmühle  
 Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392  
 Email: sales@medi-globe.de • www.medi-globe.de

# Lithotripsy Baskets – Single Use

# Lithotripsy Baskets – Single Use

# Lithotripsy Components




## Lithotripsy Basket - Braided Wire

Suitable for injectable lithotripsy spiral; tube length: 195 cm; total system length: 400 cm


Catalogue No.	Description	Sheath Outer Diameter (mm)	Suitable for Lithotripsy Spiral	Basket Width in Opened Condition (mm)	Min. Accessory Channel (mm)
GML-06-26-420	4-Wire Basket	2.6	GML-91-26-180	20	2.8
GML-06-26-425	4-Wire Basket	2.6	GML-91-26-180	25	2.8
GML-06-26-430	4-Wire Basket	2.6	GML-91-26-180	30	2.8
GML-06-26-620	6-Wire Basket	2.6	GML-91-26-180	20	2.8
GML-06-26-625	6-Wire Basket	2.6	GML-91-26-180	25	2.8
GML-06-26-630	6-Wire Basket	2.6	GML-91-26-180	30	2.8

## Lithotripsy Basket - Braided Wire

Suitable for standard lithotripsy spiral; tube length: 180 cm; total length: 400 cm



Catalogue No.	Description	Tube Ø (mm)	Suitable for Lithotripsy Spiral	Basket Width in Opened Condition (mm)	Min. Accessory Channel (mm)
GML-01-26-420	4-Wire Basket	2.6	GML-90-26-180	20	2.8
GML-01-26-425	4-Wire Basket	2.6	GML-90-26-180	25	2.8
GML-01-26-430	4-Wire Basket	2.6	GML-90-26-180	30	2.8
GML-01-26-620	6-Wire Basket	2.6	GML-90-26-180	20	2.8
GML-01-26-625	6-Wire Basket	2.6	GML-90-26-180	25	2.8
GML-01-26-630	6-Wire Basket	2.6	GML-90-26-180	30	2.8



## Lithotripsy Basket - Braided Wire

Pre-assembled in "injectable" lithotripsy spiral GML-91-26-180; total system length: 280 cm

Catalogue No.	Description	Sheath Outer Diameter (mm)	Basket Width in Opened Condition (mm)	Min. Accessory Channel (mm)
GML-11-26-420	4-Wire Basket	2.6	20	2.8
GML-11-26-425	4-Wire Basket	2.6	25	2.8
GML-11-26-430	4-Wire Basket	2.6	30	2.8
GML-11-26-620	6-Wire Basket	2.6	20	2.8
GML-11-26-625	6-Wire Basket	2.6	25	2.8
GML-11-26-630	6-Wire Basket	2.6	30	2.8

## Fitting Lithotripsy Components – Standard Lithotripsy System

Catalogue No.	Description	
GML-99-01-000	<b>StoneBuster Lithotripsy Handle</b> Reusable, user-friendly twist & pull technology, 1-push safety system	
GML-90-26-180	<b>Injectable Lithotripsy Spiral</b> Reusable, outer diameter 2.6mm, working length 180 cm	
GML-90-34-080	<b>Emergency Lithotripsy Spiral</b> Reusable, outer diameter 3.4 mm, working length 80 cm	