

SonoTip® ProControl

EUS Fine Needle Aspiration (FNA) System –

Single Use



Twist-Lock Technology

For precise tube and needle length adjustment

Laser Treated Needle Surface

For improved needle visibility

Distal Needle Stabilization

Supports tube guiding through large
working channels

Stable needle position on the target tissue

SonoTip® ProControl EUS FNA – Single Use



SonoTip® ProControl

With atraumatic (round) nitinol stylet; glass beads-blasted surface; incl. VacLok™ aspiration syringe (5–20 ml)

Catalog No.	Description	Needle Diameter (Gauge / mm)	Tube diameter (mm)	Needle length adjustment (cm)
GUS-33-18-025	EUS-Guided FNA System	25 / 0.5	1.8	0–8.5
GUS-33-27-025	EUS-Guided FNA System with tip stabilizer to accommodate large working channels (3.2 mm–3.8 mm)	25 / 0.5	2.7	0–8.5



SonoTip® ProControl

With atraumatic (round) nitinol stylet; laser engraved surface; incl. VacLok™ aspiration syringe (5–20 ml)

Catalog No.	Description	Needle Diameter (Gauge / mm)	Tube diameter (mm)	Needle length adjustment (cm)
GUS-33-18-022	EUS-Guided FNA System compatible with guide wire max. .018"	22 / 0.7	1.8	0–8.5
GUS-33-27-022	EUS-Guided FNA System with tip stabilizer to accommodate large working channels (3.2 mm–3.8 mm) compatible with guide wire max. .018"	22 / 0.7	2.7	0–8.5

SonoTip® ProControl EUS FNA – Single Use



SonoTip® ProControl

With traumatic (beveled) nitinol stylet; laser engraved surface; incl. VacLok™ aspiration syringe (5–20 ml)

Catalog No.	Description	Needle Diameter (Gauge /mm)	Tube diameter (mm)	Needle length adjustment (cm)
GUS-34-18-022	EUS-Guided FNA System compatible with guide wire max. .018"	22 / 0.7	1.8	0–8.5
GUS-34-27-022	EUS-Guided FNA System with tip stabilizer to accommodate large working channels (3.2 mm–3.8 mm) compatible with guide wire max. .018"	22 / 0.7	2.7	0–8.5



SonoTip® ProControl

With atraumatic (round) nitinol stylet; laser engraved surface; incl. VacLok™ aspiration syringe (5–20 ml)

Catalog No.	Description	Needle Diameter (Gauge /mm)	Tube diameter (mm)	Needle length adjustment (cm)
GUS-33-21-019	EUS-Guided FNA System compatible with guide wire max. .030" (recommended guide wire catalog no. GGW-08-30-400)	19 / 1.1	2.1	0–8.5



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de

Medi-Globe GmbH is a company of the Medi-Globe-Group.
For indications, contraindications, warnings and precautions, please refer to the product instructions for use.
©2024 Medi-Globe. LL-0108.EN.v11-0124