

SonoTip TopGain® – EBUS Transbronchial Needle Biopsy (TBNB) System – Single Use



Core-Tissue Sampling

For obtaining high-quality samples using an endoscopic ultrasound procedure for the primary diagnosis and staging of malignant and potentially malignant diseases

Reduction of Sample Fragmentation

Quality tissue samples through crown-cut needle tip design

Better Visibility Supports Puncture Precision

Laser engraving over the needle length

SonoTip TopGain® EBUS TBNB – Single Use



Biopsy without compromise

The new SonoTip TopGain® EBUS fine needle biopsy system was specifically developed for the high quality support of cytopathological diagnosis and to increase the quality of patient care. The innovative 3 point needle tip design with a crown cut was designed in collaboration with renowned specialists. The quality of the needle plays an important role in determining the quality of the tissue samples for the primary diagnosis and staging of malign and potentially malign diseases.



SonoTip TopGain® – compatible to Olympus (Type BF-UC190F & BF-UC180F)/Fujifilm

With atraumatic (round) Nitinol-stylet; click-lock adapter, incl. a VacLok™ aspiration syringe with a filling volume 5–20 ml

Catalogue No.	Description	Needle Material	Needle Diameter (Gauge/mm)	Needle length adjustment (cm)	Sheath Outer Diameter (mm)
GUB-45-18-022	EBUS-guided TBNB system with Twist-Lock Technology incl. separately packed adapter for Fujifilm EBUS bronchoscopes	Stainless Steel	22/0.7	0–4	1.8
GUB-46-18-022	EBUS-guided TBNB system with Twist-Lock Technology incl. separately packed adapter for Fujifilm EBUS bronchoscopes	Nitinol	22/0.7	0–4	1.8



SonoTip TopGain® – compatible to Pentax EBUS endoscopes

With atraumatic (round) Nitinol-stylet; incl. a VacLok™ aspiration syringe with a filling volume 5 – 20 ml

Catalogue No.	Description	Needle Material	Needle Diameter (Gauge / mm)	Needle length adjustment (cm)	Sheath Outer Diameter (mm)
GUB-42-18-022	EBUS-guided TBNB system with Twist-Lock Technology	Nitinol	22/0.7	0–4	1.8



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de

Medi-Globe GmbH is a company of the Medi-Globe-Group.
For indications, contraindications, warnings and precautions, please refer to the product instructions for use.
©2024 Medi-Globe. LL-0283.EN.v6-0124