

PRIMUS[®]

Advanced Single-Use Endoscopic Imaging
with PRIMUS[®] Uretero and PRIMUS[®] Cysto




- ✓ Consistently sterile
- ✓ Consistently safe
- ✓ Consistently suitable

*Your Partner in the
Endourology*

PRIMUS® Single-use flexible Ureteroscope

Single-use ureteroscope with 7.5 CH/FR outer diameter and included suction function

- ✓ Innovative
- ✓ Precise
- ✓ Versatile




- ✓ Bending lock to lock deflection of distal tip
- ✓ Ergonomic handle
- ✓ 1:1 torque
- ✓ Integrated suction function
- ✓ Rotatable sleeve for hard to reach areas and increased precision
- ✓ Fixated Tuohy Borst adapter

| Article-no. | Specifications | | | | |
|-------------|---------------------|-------------------------|----------------------|----------------------------------|---------------|
| | Outer diameter tube | Working length of shaft | Size working channel | Deflection (up/down) | Field of view |
| URS-7.5S | 7.5 CH/FR | 700 mm | 3.6 CH/FR | Standard deflection 285°/285° | 120° |
| URS-7.5S-R | | | | Reverse deflection 285°/285° | |

PRIMUS[®] Single-use flexible Cystoscope

Single-use high definition Cystoscope with 15.3 CH/FR outer diameter

- ✓ Efficient
- ✓ Reliable
- ✓ Easy handling



- ✓ 210 degree deflection in both directions
- ✓ Included Tuohy Borst adapter
- ✓ HD imaging quality
- ✓ Ergonomic handle

| Article-no. | Specifications | | | | |
|-------------|---------------------|-------------------------|----------------------|----------------------------------|---------------|
| | Outer diameter tube | Working length of shaft | Size working channel | Deflection (up/down) | Field of view |
| CYSHD-16 | 15.3 CH/FR | 400 mm | 6.6 CH/FR | Standard deflection 210°/210° | 110° |
| CYSHD-16-R | | | | Reverse deflection 210°/210° | |

PRIMUS[®] Monitor

Digital video HD-Monitor for all PRIMUS[®] endoscopic devices



- ✓ Integrated battery function
- ✓ Compatible with all PRIMUS[®] scopes
- ✓ DICOM connectivity
- ✓ Allows usage of two scopes simultaneously
- ✓ Use as stand-alone monitor or in combination with your video tower

| Article-no. | Specifications | | | | |
|-------------|-----------------|---------------|---------------------------|----------------------|--------------------|
| | Max. resolution | Display | Internal storage capacity | Digital video output | Integrated battery |
| MEDVM-D1 | 1920x1080 | 15.6" TFT LCD | 64 GB | HDMI and SDI | 4 hours |

Manufacturer

Hunan Vathin Medical Instrument Co.
1/F, Building 12,
Innovation and Entrepreneurship
Service Center,
No 9 Chuanqi West Road, Jiuhua
Economic Development Zone,
Xiangtan, Hunan, China

Exclusive Distributor

Urotech GmbH
Medi-Globe-Str. 1-5
DE-83101 Rohrdorf OT Achenmühle
Germany
Telephone: +49 (0) 8032 - 973-210
Fax: +49 (0) 8032 - 973-211
E-Mail: info@urotech.com
www.urotech.com