



Rollators & Walkers

Mobility Aids for
Everyday Use

We offer a range of rollators and walkers to suit each individual's needs. Selecting the right one to meet your requirements is an important decision that will help to increase or maintain your independence. All of our products are robust and user-friendly, promising to empower you and enable freedom on a day-to-day basis.



SAFE

With safety features such as loop lockable push down brakes, our rollators can be used with peace of mind.



FREEDOM

Designed with ease of manoeuvrability in mind.



RELIABLE

Our rollators and walkers are built from robust, high quality materials, ensuring longevity and reliability.

STEEL TRI-WHEEL WALKER

- Increase stability while walking
- Three-wheel style improves manoeuvrability without compromising support
- Loop lockable breaks with adjustable tension for safety
- Ergonomic handgrips and height-adjustable handles for comfort
- Folds easily for storage and transport
- Sturdy steel frame

Specifications

Handle height	83-93 cm
Overall depth	61 cm
Overall width	68.5 cm
Wheel diameter	19 cm
Product weight	7.7 kg

Product Code	Description
240SB	Steel Tri-Wheel Walker



ALUMINIUM LIGHTWEIGHT TRI-WHEEL WALKER

- Lightweight aluminium frame that is easy to move, making it ideal for indoor or outdoor use
- Loop lockable breaks for security and convenience
- Includes a durable bag to transport and conceal your belongings

Specifications

Handle height	83-95 cm
Overall depth	61 cm
Overall width	68.5 cm
Wheel diameter	19 cm
Product weight	7.7 kg

Product Code	Description
240L	Aluminium Lightweight Tri-Wheel



100 SERIES ROLLATOR

- Robust yet lightweight aluminium frame
- Adjustable brake tension to suit your needs
- Folds easily for storage and transport
- Ergonomic handgrips and adjustable height handles with five settings at 2.5cm intervals
- Maximum user weight of 165 kg
- Available in three different sizes to accommodate a range of heights and builds

Specifications

	Small	Medium	Large
Handle height	77.5-87.5 cm	84-94 cm	87.5-99.5 cm
Seat height	50 cm	56.5 cm	62 cm
Backrest height	75 cm	80 cm	87 cm
Wheel diameter	15 cm	20 cm	20 cm
Seat width 35 cm, Seat depth 32 cm			
Product weight	6 kg	7.2 kg	7.3 kg

Product Code	Description
102/091 566579	Small 100 Series Rollator
105/7023389	Medium 100 Series Rollator
106/091 566595	Large 100 Series Rollator



BREEZE ROLLATOR

- Alternative to a household trolley
- Made from painted aluminium tubing with four solid tyre wheels
- Includes storage bag with zip pocket, durable plastic tray and height-adjustable handles
- One-handed braking with brake lock
- Easily adjustable handle height
- Fold button for easy storage

Specifications

Handle height	83.5-94.5 cm
Overall depth	66 cm
Overall width	54 cm
Wheel diameter	17.75 cm
Tray depth	29 cm
Tray width	39 cm
Tray max load	5 kg
Bag max load	5 kg
Product weight	5.5 kg

Product Code	Description
091558360	Breeze Indoor Rollator





DELUXE LIGHTWEIGHT ROLLATOR

- Modern and compact, with enhanced manoeuvrability
- Can be folded into a 28cm standing position to fit through tight spaces, with an option to fold it further into a compact size (28cm high x 28cm wide x 39cm deep)
- Built in rest seat with canvas backrest

Specifications

Handle height	78-90 cm
Overall depth	68 cm
Seat height	53.5 cm
Backrest height	26.5 cm
Wheel diameter	20 cm
Seat width	33 cm
Seat depth	31.5 cm
Product weight	6.8 kg

Product Code	Description
091439314	Deluxe Lightweight Rollator

DAYS UPRIGHT WALKER

- Allows user to walk upright in a safer, more secure position, with adjustable height
- Cane/walking stick holder
- Padded backrest and forearm rests for comfort, can be adjusted to right length and height for user
- Made from sturdy aluminium alloy, with a maximum user weight of 150 kg
- Back of walker features two bars to assist the user in standing up

Specifications

Overall length	85 cm
Overall width	59.5 cm
Seat width	46 cm
Seat depth	28 cm
Wheel diameter	30.5 cm

Product Code	Description
7028116	Days Upright Walker



Contact Us

-  David Dunne
-  085 871 2799
-  david.dunne@fannin.eu

Fannin Ltd
Fannin House
South County Business Park
Leopardstown
Dublin 18
D18 Y0C9
T: 01 290 7000
E: customer.orders@fannin.eu
www.fannin.eu