

# Dash Carbon Power

The carbon fibre, super lightweight powerchair

The Dash Carbon Power Chair combines high strength and low weight, with a compact design that allows for ease of transportation.



- Weighs only 14.5 kg thanks to the carbon fibre frame, making the Dash Carbon Power one of the lightest power chairs on the market.
- Offset motors and detachable foot rests make the chair better for folding, for storage or putting in the car.
- Features flip up armrests and a fabric cushion for comfort.
- Built in anti-tip bars for additional safety during use.
- Optional extra battery and travel bag available.




## ENHANCED SAFETY

Built in anti-tip bars prevent chair from tipping backward.




## SUPER LIGHTWEIGHT

Weighing only 14.5 kg, making it one of the lightest power chairs on the market.



## MORE COMPACT

Offset motors and detachable foot rests allow chair to fold down to a more compact size than other power chairs.



## SPECIFICATIONS



Weight (excluding battery): 14.5 kg

Seat width	45 cm
Seat depth	40 cm
Overall width	53 cm
Overall width (closed)	24 cm
Total length	96 cm
Total height	86 cm
Rear wheel diameter	27 cm
Front wheel diameter	16.5 cm
Max user weight	130 kg

## TECHNICAL DATA

Max climable slope	10
Max speed	6 km/h
Range (medium speed on flat terrain)	> 10 km
Materials	Carbon fibre
Motor	2 x 250W brushless
Unit control	Brushless controller
Lithium battery	6.6A
Battery charger	2A



### Product Code

### Description

MSPWR0091

Dash Carbon Power Chair

## CONTACT US



David Dunne



085 871 2799



david.dunne@fannin.eu

**Fannin Ltd**  
Fannin House  
South County Business Park  
Leopardstown  
Dublin 18  
D18 Y0C9  
T: 01 290 7000  
E: customer.orders@fannin.eu  
www.fannin.eu