

# Standard E-Sac

The original Espiner Sac

## DEPLOY

- using 5mm instruments

## POSITION

- with our easy orientation tabs and black stitching to simplify identification of the Sac's mouth

## RETRIEVE


- via laparoscopic incision reducing hospital stay, cost and complication




Allowing instrument  
free extraction



One-piece tear  
resistant design with  
unique PUR coating  
for maximum strength  
and durability




Integrated  
tail allowing  
external control



# Standard E-Sac Features:



- E** Comprehensive solution covering laparoscopic procedures from node biopsy to bowel resection.
- E** Minimally invasive; reducing length of stay, financial strain and improved patient experience.
- E** Our unique coated ripstop nylon is resistant to tearing if punctured.
- E** Includes X-ray opaque marker.
- E** Polyurethane coating impermeable to biological fluids.
- E** Elongated tail providing external control throughout procedure.
- E** Easily deploy using 5mm instrument.
- E** Instrument free extraction  
– no claw forceps needed to complete.
- E** Up to a five-year shelf life.



Product Code	Trocar (mm)	Diameter (mm)	Length (mm)	Volume (ml)	Number Per Box
EMP70-INS-10	10	45	100	60	10
EMP100-INS-10	10	55	160	200	10
EMP160-INS-5	10 & 12	102	250	900	5

If you would like to know more about our original and innovative Standard E-Sac range, please call us on **+44 (0)1275 878 801** or email **[espiner@fannin.eu](mailto:espiner@fannin.eu)**  
**[www.espinermedical.com](http://www.espinermedical.com)**