

Eco Sac

The cost effective
Tissue Retrieval System

DEPLOY

- using 5mm instruments

POSITION


- with our easy orientation coloured tabs and black stitching to simplify identification of the Sac's mouth

RETRIEVE

- via laparoscopic incision reducing hospital stay, cost and complication




Safe and secure
closure with strong
monofilament
drawstring



One-piece tear
resistant design with
unique PUR coating
for maximum strength
and durability




Easy orientation
using coloured
tabs



Eco Sac Features:



- E** Comprehensive solution covering laparoscopic procedures from node biopsy to enlarged spleens.
- E** Minimally invasive; reducing length of stay, financial strain and improved patient experience.
- E** Our unique coated ripstop nylon is resistant to tearing if punctured.
- E** Includes X-ray opaque marker.
- E** Easy orientation using coloured tabs.
- E** Polyurethane coating impermeable to biological fluids.
- E** Safe closure with strong monofilament drawstring.
- E** Simplified deployment via the mouth or base.
- E** Up to a five-year shelf life.



Product Code	Trocar (mm)	Diameter (mm)	Length (mm)	Volume (ml)	Number Per Box
EMP60ECO-10	10	45	100	80	10
EMP70ECO-10	10	57	100	100	10
EMP100ECO-10	10	70	144	200	10
EMP140ECO-10	12	115	185	550	10
EMP145ECO	12	102	235	750	5
EMP170ECO	12	170	180	1200	5
EMP230ECO	12	165	340	3100	5
EMP600ECO-5	15	320	410	6000	5

To talk to Espiner about our cost effective Eco Sac range, please call us on **+44 (0)1275 878 801** or email **espiner@fannin.eu**
www.espinarmedical.com