



All-In-One ECO Video Laryngoscope

Airtraq +Plus

Blade Sizes & Geometries

- Macintosh – sizes 1 to 4
- Hyper Angulated for adults
- 90° channelled for adults, paediatrics and infants



MACINTOSH



HYPER ANGULATED



CHANNELLED

Usability

- Large 3.5" screen with rotation and image flip – enables intubation of patients in any position e.g. face-to-face. Rotation of screen can also assist intubation of large patients by shortening the handle length.
- Auto On and Auto Focus – video unit automatically sets focus for individual blade fitted
- Single use blade covers video unit handle – helps to keep unit clean
- Recording and snapshot feature – can record up to 200 videos
- VL+FoB manual light reduction

Further Usability/Teaching & Training Benefits

- Secondary screen option – enables larger screen for teaching/sharing images with colleagues
- Can be linked via physical cable to any screen with HDMI port and 640x480 image format
- Wi-Fi connectivity to secondary display or iOS/Android device
- Recording/snapshot and playback on device. Easy download via USB-C port

Sustainability

- Eco blades made from 100% recycled polycarbonate – ISCC certified (International Sustainability and Carbon Certificate)
- All packaging is made from recycled material – ISCC certified
- Rechargeable battery – purchase of replacement batteries not required



For more information:

T: 01453 874130

E: serviceuk@fannin.eu

Scan to visit our website:

