

Bacterial/Viral-HME Filters, Face Masks & accessories



- Single patient use
- Available in adult, paediatric and neonatal sizes
- Variants with elbow or catheter mounts optional
- Validated for >99.999% bacterial and viral removal
- Reduces risk of cross infection
- Latex-free

Product Code	Product Description	NHSSC Code	Units Per Box
PVC-free Face Masks			
KPFM-05	PVC-free Face Mask (large adult - size 5)	FDD3788	40
KPFM-04	PVC-free Face Mask (medium adult - size 4)	FDD3789	40
KPFM-03	PVC-free Face Mask (small adult - size 3)	FDD3790	50
KPFM-02	PVC-free Face Mask (child - size 2)	FDD3791	50
KPFM-01	PVC-free Face Mask (infant - size 1)	FDD3792	50
KPFM-00	PVC-free Face Mask (neonate - size 0)	FDD3793	50
Air Cushion Face Masks			
KACM-05	Air Cushion Face Mask (large adult - size 5)	FDD4092	25
KACM-04	Air Cushion Face Mask (medium adult - size 4)	FDD4091	25
KACM-03	Air Cushion Face Mask (small adult - size 3)	FDD4090	25
KACM-02	Air Cushion Face Mask (child - size 2)	FDD4089	50
KACM-01	Air Cushion Face Mask (infant - size 1)	FDD4088	60
KACM-00	Air Cushion Face Mask (neonate - size 0)	FDD4087	60
Bacterial/Viral-HME Filters			
KHF-10	Bacterial/Viral-HME Filter	FTCS13	50
KHFE-10	Bacterial/Viral-HME Filter + elbow	FTCS14	50
KHFX-10	Bacterial/Viral-HME Filter + expandable catheter mount 22F/straight connector 22M/15F	FTCS16	20
KHFXDE-10	Bacterial/Viral-HME Filter, expandable catheter mount 22F/straight connector 22M/15F + fixed elbow	FTCS17	20
KHFXDSE-10	Bacterial/Viral-HME Filter, expandable catheter mount 22F/straight connector 22M/15F + swivel ported elbow	FTCS18	20
KHFXE-10	Bacterial/Viral-HME Filter + expandable catheter mount 22F/fixed elbow 22M/15F	FTCS19	20
KHFXSE-10	Bacterial/Viral-HME Filter + expandable catheter mount 22F/swivel ported elbow 22M/15F	FTCS21	20
KPHF-30	Paediatric Bacterial/Viral-HME Filter	FTCS22	75
Bacterial/Viral Filters			
KF-10	Bacterial/Viral Filter	FTCS23	50
KFE-10	Bacterial/Viral Filter + elbow	FTCS25	50
KFX-10	Bacterial/Viral Filter + expandable catheter mount 22F/straight connector 22M/15F	FTCS26	20
KFXDE-10	Bacterial/Viral Filter, expandable catheter mount 22F/straight connector 22M/15F + fixed elbow	FTCS27	20
KFXDSE-10	Bacterial/Viral Filter, expandable catheter mount 22F/straight connector 22M/15F + swivel ported elbow	FTCS28	20
KFXE-10	Bacterial/Viral Filter + expandable catheter mount 22F/fixed elbow 22M/15F	FTCS29	20
KFXSE-10	Bacterial/Viral Filter + expandable catheter mount 22F/swivel ported elbow 22M/15F	FTCS31	20
Catheter mounts and elbow pieces			
KX-10	Expandable catheter mount 22F/straight connector 22M/15F	FDB719	50
KXE-10	Expandable catheter mount 22F/fixed elbow 22M/15F	FDB547	50
KXSE-10	Expandable catheter mount 22F/swivel ported elbow 22M/15F	FDB721	50
KE-10	Fixed elbow piece	FDC779	100
KSE-10	Swivel ported elbow piece	FDB722	100
Airway Component Packs – Bacterial/Viral-HME Filters			
KTHFXDE-05	Air Cushion Face Mask (large adult - size 5), Bacterial/Viral-HME Filter, expandable catheter mount 22F/straight connector 22M/15F + fixed elbow	FDD3771	20
KTHFXDE-04	Air Cushion Face Mask (medium adult - size 4), Bacterial/Viral-HME Filter, expandable catheter mount 22F/straight connector 22M/15F + fixed elbow	FDD3772	20
KTHFXDE-03	Air Cushion Face Mask (small adult - size 3), Bacterial/Viral-HME Filter, expandable catheter mount 22F/straight connector 22M/15F + fixed elbow	FDD3773	20
KTHFXE-05	Air Cushion Face Mask (large adult - size 5), Bacterial/Viral-HME Filter, expandable catheter mount 22F/fixed elbow 22M/15F	FDD3774	20
KTHFXE-04	Air Cushion Face Mask (medium adult - size 4), Bacterial/Viral-HME Filter, expandable catheter mount 22F/fixed elbow 22M/15F	FDD3775	20
KTHFXE-03	Air Cushion Face Mask (small adult - size 3), Bacterial/Viral-HME Filter, expandable catheter mount 22F/fixed elbow 22M/15F	FDD3776	20