

Double Lumen Silicone Laryngeal Masks

- Reusable
- Reinforced
- Medical grade silicone
- Validated for up to 40 uses
- DEHP-free
- Latex free





Ordering details

Product Code	Size	Patient Weight	Cuff Volume	Units per Box	NHSSC Code
FF30810	1	< 5 kg	< 4ml	5	FDD85143
FF30815	1.5	5-10 kg	< 7ml	5	FDD85153
FF30820	2	10-20 kg	< 10ml	5	FDD85141
FF30825	2.5	20-30 kg	< 14ml	5	FDD85137
FF30830	3	30-50 kg	< 20ml	5	FDD85167
FF30840	4	50-70 kg	< 30ml	5	FDD85176
FF30850	5	70-100 kg	< 40ml	5	FDD85182



Fannin Ltd
 Fernhouse - Suite 3
 Daniels Industrial Estate
 Bath Road, Gloucestershire
 GL5 3TJ, United Kingdom
T: +44 (0) 1453 874 130
E: serviceuk@fannin.eu