


# Tracheal tube introducers

- Tracheal tube introducers may be used to facilitate endotracheal intubation
- Made of medical grade PE, latex-free, flexible
- Distal rounded smooth tip allows introducer to be inserted, without damaging soft tissues
- Memory retention allows shaping to clinical needs
- Clearly graduated markings provide depth indication, to ensure correct placement
- Different sizes available, suitable for paediatric and adult use
- Single use




Straight tip




Curved tip



Vented curved











Vented straight



## Difficult intubations

- Tracheal tube introducers may be used during difficult intubations

Visualised structures	Grade 1	Grade 2	Grade 3	Grade 4
Pharyngees (Mallanpati)				
Laryngees (Cormack & Lehane)				

If greater than Stage 2, Mallanpati risk of DI due to bad visualisation of vocal chords

## Non-vented



Product code	Product description	NHSSC code	Units per box
FFI0006	Tracheal tube introducer, curved tip, 6Fr, 53.5cm	FSM85378	10
FFI0008	Tracheal tube introducer, curved tip, 8Fr, 53.5cm	FSM85384	10
FFI0014	Tracheal tube introducer, curved tip, 14Fr, 70cm	FSM85381	10
FFI0106	Tracheal tube introducer, straight, 6Fr, 53.5cm	FSM85376	10
FFI0108	Tracheal tube introducer, straight, 8Fr, 53.5cm	FSM85388	10
FFI0114	Tracheal tube introducer, straight, 14Fr, 70cm	FSM85389	10



## Coiled

- Innovative, compact packaging, designed specifically for use in pre-hospital



Product code	Product description	NHSSC code	Units per box
FFI0310	Tracheal tube introducer, coiled, curved tip, 10Fr, 70cm	FSM85377	10
FFI0315	Tracheal tube introducer, coiled, curved tip, 15Fr, 70cm	FSM85387	10



## Vented



Product code	Product description	NHSSC code	Units per box
FFI0208	Tracheal tube introducer, vented, curved tip, 8Fr, 53.5cm	FSM85390	5
FFI0210	Tracheal tube introducer, vented, curved tip, 10Fr, 70cm	FSM85379	5
FFI0214	Tracheal tube introducer, vented, curved tip, 14Fr, 70cm	FSM85382	5
FFI0208U	Tracheal tube introducer, vented, curved tip, 8Fr, 53.5cm, without 15mm connector	—	5
FFI0210U	Tracheal tube introducer, vented, curved tip, 10Fr, 70cm, without 15mm connector	—	5
FFI0214U	Tracheal tube introducer, vented, curved tip, 14Fr, 70cm, without 15mm connector	—	5

# Exchanger

- Used to facilitate uncomplicated, endotracheal tube exchanges
- Made of medical grade PE, latex-free, flexible
- Atraumatic tip
- Distal side port helps to ensure sufficient air flow
- Clearly graduated marking provide depth indication, to ensure correct placement
- Removable connector permits ventilation if necessary during exchange procedure



Product code	Product description	NHSSC code	Units per box
FFI0708	Tracheal tube introducer, exchanger, 8Fr, 53.5cm	—	5
FFI0710	Tracheal tube introducer, exchanger, 10Fr, 70cm	—	5
FFI0714	Tracheal tube introducer, exchanger, 14Fr, 70cm	—	5

