

# NacSax®

## VOMMAX®

Simple to use and effective Vomit Bag



9310-002

**Vommax®** is the ideal replacement for the traditional disposable bowl and kidney dish - Ergonomically designed with a rigid plastic mouthpiece fitting perfectly over the patients mouth ensuring that all vomit is completely captured, thus, eliminating accidents.

**Vommax®** has a simple twist and seal device for immediate hygienic odourless disposal

**Vommax®** is manufactured using a translucent material, allowing any blood to be easily identified and indicators allowing approximate volume assessment.

**Vommax®** is supplied with or without a solidifying gel agent.

**Vommax®** is used in many areas including;

- Hospitals
- Ambulances
- A & E Departments
- Doctor's Surgeries
- Any Transport situation where there could be a possibility of travel sickness

# VacSax®

## VOMMAX®



GRIP



TWIST



LOCK



DISPOSE

9310-003

Part No.	Description
9310-002	<b>Vommax®</b> Pack of 50
9310-003	<b>Vommax®</b> with solidifying gel agent. Pack of 50
9310-004	Perspex Holder to accept Pack of 50 <b>Vommax®</b>

### VacSax Limited

Western Wood Way, Langage Science Park,  
Plymouth, Devon PL7 5BG, U.K.

**Tel:** +44 (0)1752 337000 **Fax:** +44 (0)1752 337600

**Email:** [sales@vac.sax.com](mailto:sales@vac.sax.com)

**[www.vacsax.com](http://www.vacsax.com)**

Our **commitment** is to **quality**