



INSUFFLATION TUBING KITS & PARTICULATE FILTERS

- Simple and easy way to keep your patient and valuable equipment safe.
- Feel confident that the gas you are delivering to your patient's abdomen is free from debris and contaminants.
- Eliminate cross contamination between patients and prevent the backflow of blood and body fluids entering the insufflator



INSUFFLATION TUBING KITS & PARTICULATE FILTERS

- **Cost effective**
- **Sterile** – Single use – Easy peel pouch
- **Bi Directional Bacterial/Viral Filter** – 99.9999% efficiency at 0.027 micron
- **Excellent Performance** – High flow rate – Low Resistance – Enables rapid inflation
- **Flexible 3m tubing**
- **Luer Lock Connector**, Wolf Connector or push fit connector
- **Latex Free**



IFT0001AD



IFT0002AD



IFT0015AD

Fannin - giving you choice quality and savings

Skintact insufflation kits are designed to meet the needs of a busy theatre environment. Reduce preparation time for laparoscopic procedures by allowing to provide you with a combination of products in one convenient pouch.

Each item is sterilized separately giving you the flexibility of using all of the components or just the ones you choose for that particular case. The unused items remain sterile and can be retained for future use.

Kits are customisable to suit your particular needs.

Ordering details

Skintact Insufflation Filter Kits

Code	Description	UOI	NHS Supply Chain code
IFT0003	Insufflation Wolf IFT0002+ Easi-Drape 59143 + Easi-C Anti-Fog	10	
IFT0005	Insufflation Luer Lock IFT0001+ Needle Counter NC1010 + Easi-Drape 59143	10	
IFT0006	Insufflation Luer Lock IFT0001+ Easi-Drape 59143	20	FSW349
IFT0007	Insufflation Luer Lock IFT0001 + Easi-C Anti-Fog Easi-Drape 59143	10	FSW350
IFT0008	Insufflation Wolf IFT0002+ Easi-Drape 59143	10	
IFT0010	Insufflation Luer Lock IFT 0001+ Easi-Drape 59144 + Easi-C Anti -fog	20	
IFT0012	Insufflation Luer Lock IFT0001+ Easi-Drape 59143 + Easi-C Anti-Fog + Passive Smoke Filter PSF 0004	10	
IFT0013	Insufflation Luer Lock IFT0001+ Passive Smoke Filter PSF 0004	10	
IFT0014	Insufflation Luer Lock IFT0001 + Easi-C Anti-Fog FC00013	10	
IFT0017	Insufflation Push fit IFT0015 + 59143 Easi-Drape	10	
IFT0018	Insufflation Push fit IFT0015 + 59143 Easi-Drape + FC0013 Easi-C	10	
IFT0024	Lap Hernia Pack - IFT0001 (Insufflation Tubing), 59143 (Camera Drape), FC00013 (Anti-Fog), NC0001 (Sharps Pad), 20ml Syringe & 20g Needle	10	
IFT0025	Insufflation Push IFT0025 + Easi-C Anti-Fog	10	

Skintact Insufflation Tubing /Filter

IFT0001AD	Insufflation Filter with male Luer connectors both ends	25	FSW346
IFT0002AD	Insufflation Filter with 15mm Wolf Connector	25	FSW347
IFT0015AD	Insufflation Filter with Push fit connector	25	FSW348