


# Laparoscopic

## Liquid Scope Warmer

**Mediflex®**  
SURGICAL PRODUCTS



### Liquid Scope Warmer

A simple, pre-surgical Scope Warmer which greatly reduces fogging during procedures. The self-contained liquid pack provides rapid warming without the waste of time and cost associated with alternative anti-fogging solutions.

- Simple activation to start warming action
- Quick process, warms scopes in under 1 minute
- Provides lens preservation – no hard surface or chemicals
- Ideal for robotic cases
- Accommodates all scope sizes
- Eliminates 'hot water' approach
- Cost-effective, inexpensive
- Minimal storage space required



### ORDERING INFORMATION

Ref. No.	Description
74005-20	Liquid Scope Warmer, Sterile, Box of 20 <i>Dimensions (approx.): 7" x 2" (18cm x 5cm)</i>

**For more information, please contact us at:**

**Fannin UK Ltd**  
Fernhouse - Suite 3  
Daniels Industrial Estate  
Bath Road, Gloucestershire  
GL5 3TJ, United Kingdom  
**T:** +44 (0) 1453 874 130  
**E:** salesll@fannin.eu

**Fannin Ltd.**  
Fannin House  
South County Business Park  
Leopards Town  
Dublin 18  
**T:** +353 (0)1 290 7000  
**E:** info@fannin.eu

[www.fannin.eu](http://www.fannin.eu)