



SofTip™ Microlaryngeal Tubes

SofTip™ technology and AccuCuff™ cuff pressure indicators provide dual protection for patient airways




AccuCuff™ Cuff Pressure Indicator

Permits real time monitoring of the cuff pressure, avoiding excessive or insufficient air inflation




SofTip™ Technology

Allows for faster and safer intubations, reducing the trauma to patient airways



High Flexibility

Made of flexible and soft material



Smaller Diameter & Longer Length

Allows for better access to the laryngeal area, enables surgeons to have a clear view and provide more space for the surgical team to work



Features

- Specifically designed for laryngeal procedures and microlaryngeal surgery
- Serves as an alternative when larger size ETTs prove difficult to intubate
- Latex and DEHP free
- 5-year shelf life

Ordering Details

SofTip™ Microlaryngeal Tube, Oral/Nasal with AccuCuff™, HVLP Cuff

Size	ID. mm \pm 0.15	OD. mm \pm 0.15	Cuff mm \pm 0.15	Length mm	Code	Units per Box
4.0mm	4.0	5.3	28	360	ETT-YCP32-40	10
5.0mm	5.0	6.7	28	360	ETT-YCP32-50	10
6.0mm	6.0	8.0	28	360	ETT-YCP32-60	10



NHSSC product codes for ordering

Scan the QR code:



Fannin Ltd
 Fernhouse - Suite 3
 Daniels Industrial Estate
 Bath Road, Gloucestershire
 GL5 3TJ, United Kingdom
 T: +44 (0) 1453 874 130
 E: serviceuk@fannin.eu