



SofTip™ Endotracheal Tubes

SofTip™ technology can make the entire intubation process faster and safer



- Atraumatic, soft and flexible
- Extensive range of standard, reinforced and preformed cuffed tubes available
- High volume, low pressure cuff
- Clear connector – EGBAT compliant
- Murphy eye
- Single use
- Latex free and DEHP free
- 5-year shelf life



SofTip™ technology



High quality ISO connectors

Ordering Details

Standard

Size	Code	Size	Code	Units per box
2.5mm	ETT-P22-25*	6.5mm	ETT-P22-65	10
3.0mm	ETT-P22-30	7.0mm	ETT-P22-70	10
3.5mm	ETT-P22-35	7.5mm	ETT-P22-75	10
4.0mm	ETT-P22-40	8.0mm	ETT-P22-80	10
4.5mm	ETT-P22-45	8.5mm	ETT-P22-85	10
5.0mm	ETT-P22-50	9.0mm	ETT-P22-90	10
5.5mm	ETT-P22-50	9.5mm	ETT-P22-95	10
6.0mm	ETT-P22-60	10.0mm	ETT-P22-100	10

* SofTip™ technology not available for this size tip

South Facing

Size	Code	Size	Code	Units per box
2.5mm	-	6.5mm	ETT-K22-65	10
3.0mm	ETT-K22-30	7.0mm	ETT-K22-70	10
3.5mm	ETT-K22-35	7.5mm	ETT-K22-75	10
4.0mm	ETT-K22-40	8.0mm	ETT-K22-80	10
4.5mm	ETT-K22-45	8.5mm	ETT-K22-85	10
5.0mm	ETT-K22-50	9.0mm	ETT-K22-90	10
5.5mm	ETT-K22-50	9.5mm	ETT-K22-95	10
6.0mm	ETT-K22-60	10.0mm	ETT-K22-100	10

Reinforced

Size	Code	Size	Code	Units per box
2.5mm	-	6.5mm	ETT-J22-65	10
3.0mm	ETT-J22-30	7.0mm	ETT-J22-70	10
3.5mm	ETT-J22-35	7.5mm	ETT-J22-75	10
4.0mm	ETT-J22-40	8.0mm	ETT-J22-80	10
4.5mm	ETT-J22-45	8.5mm	ETT-J22-85	10
5.0mm	ETT-J22-50	9.0mm	ETT-J22-90	10
5.5mm	ETT-J22-50	9.5mm	ETT-J22-95	10
6.0mm	ETT-J22-60	10.0mm	ETT-J22-100	10

North Facing

Size	Code	Size	Code	Units per box
2.5mm	-	6.5mm	ETT-B22-65	10
3.0mm	ETT-B22-30	7.0mm	ETT-B22-70	10
3.5mm	ETT-B22-35	7.5mm	ETT-B22-75	10
4.0mm	ETT-B22-40	8.0mm	ETT-B22-80	10
4.5mm	ETT-B22-45	8.5mm	ETT-B22-85	10
5.0mm	ETT-B22-50	9.0mm	ETT-B22-90	10
5.5mm	ETT-B22-50	9.5mm	ETT-B22-95	10
6.0mm	ETT-B22-60	10.0mm	ETT-B22-100	10

NHSC product codes for ordering

Scan the QR code:



Fannin Ltd
Fernhouse - Suite 3
Daniels Industrial Estate
Bath Road, Gloucestershire
GL5 3TJ, United Kingdom
T: +44 (0) 1453 874 130
E: serviceuk@fannin.eu
fannin.eu

