


ArtCan Arterial Safety Cannula

- Suitable for arterial pressure monitoring and serial arterial blood gas determinations
- Sharp atraumatic needle bevel for minimal patient discomfort
- Specially designed ON/OFF switch prevents backflow and exposure of blood to healthcare professionals, and provides patient safety against air embolism and contamination
- Needle safety guard automatically covers the sharp bevel, after withdrawal of needle from the hub, minimising the risk of needle stick injury
- Available with PUR catheter, 45mm
- Sizes: 20G & 22G
- Latex-free, non-DEHP and PVC-free



ArtCan Arterial Safety Cannula



ON/OFF switch



Needle safety guard

Ordering information

Product code	Product description	Size	Flow rate (mL/min)	Units per box	NHSSC code
13425	ArtCan Arterial Safety Cannula, PUR catheter	20G	49	25	FSP85890
13429	ArtCan Arterial Safety Cannula, PUR catheter	22G	28	25	FSQ3473

Sterilisation Method: Ethylene Oxide sterilised