


Autofusion Enhance IV Sets




- Autofusion Enhance IV Sets with air vents and 15µm filters
- Hydrophilic membrane to prevent air entrainment if IV bag empties
- Available with 3-way connector and extension tubing
- Needle-free valves as optional extras


Autofusion Enhance IV Sets




1029DFAFE
Autofusion Enhance IV Set with air vent, 180cm tubing, rotating luer connector




1029DF2AFE
Autofusion Enhance IV Set with air vent, 240cm tubing, rotating luer connector




1029DFVAFE
Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve




1029DFXAFE
Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector, 120cm extension tubing, rotating luer connector




1029DFXVAFE
Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve, 3-way connector, 120cm extension tubing, rotating luer connector




1029DFXMYAFE
Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector with needle-free valve, 120cm extension tubing, rotating luer connector



1029DFXVMYAFE
Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve, 3-way connector, with needle-free valve, 120cm extension tubing, rotating luer connector



1048BMYAFE
Autofusion Enhance IV Set with air vent, 210cm tubing, needle free access port, back check valve, rotating luer connector



1048B2MYAFE
Autofusion Enhance IV Set with air vent, 260cm tubing, two needle free access ports, back check valve, rotating luer connector

Product code	Product description	Units per box	NHS Supply Chain code
1029DFAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, rotating luer connector	50	FSB2113
1029DF2AFE	Autofusion Enhance IV Set with air vent, 240cm tubing, rotating luer connector	50	FSB85130
1029DFVAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve	50	FSB2114
1029DFXAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector, 120cm extension tubing, rotating luer connector	50	FSB2115
1029DFXVAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve, 3-way connector, 120cm extension tubing, rotating luer connector	50	FSB2116
1029DFXMYAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector, with needle-free valve, 120cm extension tubing, rotating luer connector	50	FSB2117
1029DFXVMYAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve, 3-way connector, with needle-free valve, 120cm extension tubing, rotating luer connector	50	FBO85003
1048BMYAFE	Autofusion Enhance IV Set with air vent, 210cm tubing, needle free access port, back check valve, rotating luer connector	50	FBO85001
1048B2MYAFE	Autofusion Enhance IV Set with air vent, 260cm tubing, two needle free access ports, back check valve, rotating luer connector	50	FBO85004