

Integral Silicone Laryngeal Masks

Standard & Reinforced

- Single use
- 100% medical grade silicone
- Colour-coded inflation balloons for size identification
- Available in all sizes
- PVC-free
- Latex-free



High quality ISO
connections



Highly flexible
silicone tube



Laryngeal mask
introducer



Silicone cuff enables
atraumatic insertion
& removal

Ordering details

Standard Integral Silicone Laryngeal Masks - single use

Product Code	Size	Units per Box	Recommended Patient Weight	Inflation Valve Colour	NHSSC Code
LM1004S2	1	10	Neonates/infants <5 kg	 Purple	FDD1000
LM1507S2	1.5	10	Infants 5-10 kg	 White	FDD1001
LM2010S2	2	10	Infants/children 10-20 kg	 Red	FDD1002
LM2514S2	2.5	10	Children 20-30 kg	 Orange	FDD1003
LM3020S2	3	10	Children 30-50 kg	 Green	FDD1004
LM4030S2	4	10	Adults 50-70 kg	 Yellow	FDD1005
LM5040S2	5	10	Adults 70-100 kg	 Pink	FDD1006

Reinforced Integral Silicone Laryngeal Masks - single use

Product Code	Size	Units per Box	Recommended Patient Weight	Inflation Valve Colour	NHSSC Code
LM1004S2R	1	10	Neonates/infants <5 kg	 Purple	FDD1007
LM1507S2R	1.5	10	Infants 5-10 kg	 White	FDD1008
LM2010S2R	2	10	Infants/children 10-20 kg	 Red	FDD1009
LM2514S2R	2.5	10	Children 20-30 kg	 Orange	FDD1010
LM3020S2R	3	10	Children 30-50 kg	 Green	FDD1011
LM4030S2R	4	10	Adults 50-70 kg	 Yellow	FDD1012
LM5040S2R	5	10	Adults 70-100 kg	 Pink	FDD1013

Laryngeal Mask Introducer - single use

Product Code	Units per Box	Recommended Patient Weight	NHSSC Code
INT304	10	Medium 3 & 4	FDD1720