

Transoesophageal Echo Transducer Cover Kit

Sterile Ultrasound Transducer Cover and Kits for Cardiology TOE/TEE Transducers

- Sterile single use, with high quality and durable Polyurethane
- Specifically designed for TOE transducers: Tip: 16mm over 400 mm in length, tapered to 31 mm over 700 mm in length | **Total length: 1100 mm**
- Unique 'Easi' applicator - Quick and efficient insertion of ultrasonic gel and smooth release of cover - Applicator fits securely onto most control handles
- Extended 10ml Syringe pre-loaded with Ultrasonic Transmission Gel
- Bite guard included
- Latex-free
- 3 year shelf life

Ordering Details

Product Code	Product Description	Dimension	Box Size
Adult 1250-02	Cardiology TOE/TEE covers Protection™ TEE kit sterile – Square Tip	15/31 x 1100mm	24
Adult 1251-01	Cardiology TOE/TEE covers Protection™ TEE kit sterile – Round Tip	15/31 x 1100mm	24
Paediatric 1260-02	Cardiology TOE/TEE covers Protection™ TEE kit sterile – Round Tip	11/31 x 1000mm	24
Neonatal 1280-02	Cardiology TOE/TEE covers Protection™ TEE kit sterile – Round Tip	11/31 x 1000mm	24

