

Transoesophageal Echo Transducer Cover Kit

Sterile Ultrasound Transducer Cover and Kits for Cardiology TOE/TEE Transducers

- Sterile single use, with high quality and durable Polyurethane
- Specifically designed for TOE transducers: Tip: 16mm over 400 mm in length, tapered to 31 mm over 700 mm in length | **Total length: 1100 mm**
- Unique 'Easi' applicator - Quick and efficient insertion of ultrasonic gel and smooth release of cover - Applicator fits securely onto most control handles
- Extended 10ml Syringe pre-loaded with Ultrasonic Transmission Gel
- Bite guard included
- Latex-free
- 3 year shelf life

Ordering Details

| Product Code | Product Description | Dimension | Box Size |
|--------------------|---|----------------|----------|
| Adult 1250-02 | Cardiology TOE/TEE covers Probtection™ TEE kit sterile – Square Tip | 15/31 x 1100mm | 24 |
| Adult 1251-01 | Cardiology TOE/TEE covers Probtection™ TEE kit sterile – Round Tip | 15/31 x 1100mm | 24 |
| Paediatric 1260-02 | Cardiology TOE/TEE covers Probtection™ TEE kit sterile – Round Tip | 11/31 x 1000mm | 24 |
| Neonatal 1280-02 | Cardiology TOE/TEE covers Probtection™ TEE kit sterile – Round Tip | 11/31 x 1000mm | 24 |

