

Ultrasound Probe Covers

Disposable ultrasound probe covers designed to help maintain a sterile field during procedures

- Sterile, single-use probe covers, manufactured from high quality and durable polyurethane
- Available in many shapes and sizes, either as a separate cover or as a complete set including ultrasound transmission gel
- Helps protect and may extend equipment and instrument life. Provides an effective contagion barrier for your equipment.

Ordering Details

Product Code	Product Description	Dimension	Box Size	NHSS Code
59156	Disposable Non Latex Polyethylene Ultrasound Probe Cover + 2 Fixation Bands	150 mm x 0.85 m	100	FBEI0146
59157	Disposable Non Latex Polyethylene Ultrasound Probe Cover + 2 Fixation Bands 150 mm x 0.85 m (59156) + Sterile Gel	150 mm x 0.85 m	50	FBEI0147
59158	Disposable Non Latex Universal Polyethylene Ultrasound Probe Cover With Ring Applicator + 2 Fixation Bands	150 mm x 2.5 m	50	FBEI2389
59159	Disposable Non Latex Universal Polyethylene Ultrasound Probe Cover With Ring Applicator + 2 Fixation Bands + Sterile Gel	150 mm x 2.5 m	50	Available Direct
59166	Disposable Non Latex Universal Ring Camera Drape 150 mm x 1.275 m + (1236-02) Polyurethane Laparoscopic Ultrasound Cover 22 mm x 540 mm & Twist Lock	150 mm x 1.275 m	100	Available Direct
59167	Disposable Non Latex Polyethylene Ultrasound Probe Cover + 2 Fixation Bands 150 mm x 1.275 m (59166) + Sterile 20 g Sachet of Ultrasound Gel	150 mm x 1.275 m	50	Available Direct
1205-01	Non Latex Non-Sterile Polyurethane Probe Cover & Twist Lock	30 mm x 400 mm	50	Available Direct



Fannin Ltd

Fernhouse - Suite 3 Daniels
Industrial Estate Bath Road,
Gloucestershire GL5 3TJ,
United Kingdom
T: +44 (0) 1453 874 130
E: serviceuk@fannin.eu
fannin.eu

Scan to view our
online catalogue:

