

# Ultrasound Transducer Covers

## Single-Use Transducer Covers for Gynaecology and Urology

- Single-use with high quality and durable polyurethane
- Latex and latex-free options available
- Optimal fit for all of the trans-vaginal and trans-rectal transducers currently in use
- Individually packed
- Helps protect patients and healthcare workers against cross-contamination
- Specifically designed for trans-vaginal transducers: Tip: 30 mm over 300 mm in length
- Available across TOE, transvaginal and ultrasound applications and in sterile and non-sterile configurations

### Ordering Details

| Product Code | Product Description   | Dimension      | Box Size |
|--------------|---|----------------|----------|
| 1000LL       | Non Latex non-sterile rolled condom style polyisoprene probe cover      | 33 mm x 195 mm | 100      |
| 1005         | Non Latex non-sterile rolled condom style polyisoprene probe cover      | 26 mm x 200 mm | 100      |
| 1006LL       | Non Latex non-sterile rolled condom style polyisoprene probe cover      | 30 mm x 300 mm | 100      |
| 1100LL       | Latex non-sterile rolled condom style probe cover Gynaecology / Urology | 33 mm x 195 mm | 144      |
| 1104         | Latex non-sterile probe cover Gynaecology / Urology                     | 26 mm x 300 mm | 100      |
| 1106         | Latex non-sterile rolled condom style probe cover                       | 30 mm x 300 mm | 100      |

