

Transoesophageal Echo Transducer Cover Kit

Sterile, Single-Use Ultrasound Transducer Cover and Kits for Cardiology TOE/TEE Transducers

- Sterile, single-use TOE/TEE probe covers manufactured from high quality, durable polyurethane
- Tip: 16 mm over 400 mm in length, tapered to 31 mm over 700 mm in length | **Total length: 1100 mm**
- Unique 'Easi' applicator - quick and efficient insertion of ultrasound transmission gel and smooth release of cover. Applicator fits securely onto most control handles
- Extended 10 ml syringe pre-loaded with ultrasound transmission gel
- Bite guard included
- Latex-free
- 3-year shelf life
- Helps protect patients and healthcare workers against cross-contamination

Ordering Details

Product Code	Product Description	Dimension	Box Size
Adult 1250-02	Cardiology TOE/TEE covers Protection™ TEE kit sterile – Square Tip	15/31 mm x 1100 mm	24
Adult 1251-01	Cardiology TOE/TEE covers Protection™ TEE kit sterile – Round Tip	15/31 mm x 1100 mm	24
Paediatric 1260-02	Cardiology TOE/TEE covers Protection™ TEE kit sterile – Round Tip	11/31 mm x 1000 mm	24
Neonatal 1280-02	Cardiology TOE/TEE covers Protection™ TEE kit sterile – Round Tip	11/31 mm x 1000 mm	24

