

Non-Sterile Apparel

Fannin UK offers a comprehensive range of disposable personal protective apparel designed for use across operating theatres, clinical environments and healthcare settings.

Disposable Warm-up Jackets

Disposable warm-up jackets provide lightweight warmth and comfort in cooler clinical environments. Manufactured from breathable SMS material with knitted cuffs and snap fastenings, they help maintain comfort while supporting infection prevention.

Facial Protection

Our range of surgical and procedural masks provides dependable protection without compromising comfort or breathability, helping support staff throughout routine and extended procedures.

Surgical Gowns & Scrub Suits

Available in AAMI Levels 1 - 4, our surgical gowns offer appropriate fluid and microbial barrier protection for procedures ranging from minimal-risk examinations to high-fluid surgical interventions. Our disposable scrub suits offer comfort, breathability and freedom of movement while supporting high standards of hygiene. Each set includes a V-neck scrub top and matching trousers.

Shoe Covers

Disposable shoe covers help maintain clean, controlled clinical environments by reducing contamination. Available in polypropylene and polyethylene options, they provide a secure, comfortable fit for everyday healthcare use.

Surgical Headwear

Our range of disposable surgical caps and mob caps is designed to help maintain sterile environments while providing a comfortable, secure fit for healthcare professionals.



Fannin Ltd

Fernhouse - Suite 3 Daniels
Industrial Estate Bath Road,
Gloucestershire GL5 3TJ,
United Kingdom
T: +44 (0) 1453 874 130
E: serviceuk@fannin.eu
fannin.eu

Scan to view our online
catalogue:



Ordering Information

Warm Up Jackets

Product code	Product Description	Unit of Issue	Size	Ordering Information
WG0084LN-IH	Disposable Warm Up Jacket	50	Large	Available Direct via Fannin UK
WG0084MN-IH	Disposable Warm Up Jacket	50	Medium	Available Direct via Fannin UK
WG0084XLN-IH	Disposable Warm Up Jacket	50	Extra Large	Available Direct via Fannin UK

Shoe Covers

Product code	Product Description	Unit of Issue	Size	Ordering Information
4850-IH	Disposable Shoe Covers	1500	Large	Available Direct via Fannin UK
4860-IH	Disposable Shoe Covers	1500	Regular	Available Direct via Fannin UK

Surgical Headwear

Product code	Product Description	Unit of Issue	Size	Ordering Information
3274-IH	Surgical Cap Bouffant (Blue)	1500	Large	Available Direct via Fannin UK
3244-IH	Surgical Cap Bouffant (White)	2000	Regular	Available Direct via Fannin UK
4359-IH	Surgical Cap 40G (Blue)	200	Universal	Available Direct via Fannin UK

Scrubs

Product code	Product Description	Unit of Issue	Size	Ordering Information
90506-IH	Scrub Suit, Blue	50	Small	Available Direct via Fannin UK
90507-IH	Scrub Suit, Blue	50	Medium	Available Direct via Fannin UK
90508-IH	Scrub Suit, Blue	50	Large	Available Direct via Fannin UK
90509-IH	Scrub Suit, Blue	50	Extra Large	Available Direct via Fannin UK
90510-IH	Scrub Suit, Blue	50	Extra Extra Large	Available Direct via Fannin UK



Contact your Fannin UK representative for more information or to place an order.

Fannin Ltd

Fernhouse - Suite 3 Daniels
Industrial Estate Bath Road,
Gloucestershire GL5 3TJ,
United Kingdom

T: +44 (0) 1453 874 130

E: serviceuk@fannin.eu

fannin.eu

Scan to view our
online catalogue:

